

Patient Transfer

3 12

EL

-



Improving life.

Discover why the EZ-Tranz Gold Patient Transfer Board is the best solution for lateral patient transfer



ΕZ

Easy to use because of the instruction diagram. Wall rack available to keep it nearby.

SAFE

Safer for staff because there is no need to lift patients to transfer. Only 3 lbs. EZ to lift and position, reducing fatigue and injury to staff. MRI safe!

COMFORTABLE

Patients roll quietly and smoothly on the solid foam core which is much more comfortable than other devices.

CLEAN

Smooth Antimicrobial Plastic handles coupled with an Ultra Fresh antimicrobial cover make this board a great choice in the fight against HAIs. Can be cleaned with most antibacterial wipes found in healthcare facilities.

MCM190 Series

EZ-Tranz Gold Patient Transfer Board

Gold Antimicrobial cover around a soft durable core, make this EZ-Tranz Board an elite transfer aid. Conveyor type action allows for smooth lateral transfer of patients. Extremely lightweight and easy to use make this the first choice for medical professionals and patients.

- > Easy to clean
- Antimicrobial
- > Extremely lightweight

Replacement covers fit corresponding length of transfer board. Not guaranteed to fit other brands of transfer boards.

MCM195 30"L Replacement Cover MCM196 48"L Replacement Cover

MCM123

MCM125 MCM130



MCM120 Series Patient Transfer Roller Boards

Roller style patient transfer boards easily transfer patients from bed to gurney or OR table. Constructed of aluminium for easy to clean and minimal weight reasons. Standard style patient transfer board. 15" width and 3 lengths to meet your needs.

Replacement covers fit corresponding length of roller board. Not guaranteed to fit other brands of roller boards.

15"W x 25"L

15"W x 30"L

15"W x 67"L

| MCM124 | 25"L Replacement Cover |
|--------|---------------------------|
| MCM126 | 30"L Replacement Cover |
| MCM131 | 67″L Replacement Cover |



MCM120 Series

Slider Style Patient Transfer Boards

Constructed of HDPE. 3/16" thick. Non-anti-static. Generous handle for easy transfers. Easily transfer patients from bed to gurney or OR table. Smooth easy to clean surface. Available in two standard sizes.

 MCM120
 22"W x 33"L x 3/16"H

 MCM122
 22"W x 72"L x 3/16"H

MCM190 15"W x 30"L

MCM191 15"W x 48"L



MCM190 Series

. . . .

Wall Rack for EZ-Tranz Gold Patient Transfer Boards

MCM199 E

MCM132

Wall Rack for EZ-Tranz Gold Patient Transfer Board

Constructed of stainless steel material for ease of cleaning and prevention of infection spread. Fits all sizes of EZ-Tranz Gold transfer boards. Keeps boards clean and out of the way, yet in the spot you need them. Each rack holds one board. Mounting hardware not included due to varying substrates.



Constructed of stainless steel material for ease of cleaning and prevention of infection spread. Fits all types of roller boards. Keeps boards clean and out of the way, yet in the spot you need them. Each rack holds one board. Mounting hardware not included due to varying substrates.



MIDCENTRAL® MEDICAL

MCM120 Series

Wall Rack for Slider Style Patient Transfer Boards

Constructed of stainless steel material for ease of cleaning and prevention of infection spread. Hook holds small or large size patient transfer slider boards. Keeps boards clean and out of the way, yet in the spot you need them. Each rack holds one board. Mounting hardware not included due to varying substrates.

Wall Rack for MCM120 Slider Style Patient Transfer Boards

Wall Rack for Patient Transfer Roller Boards

